

Environmental Health & Safety Policy Manual									
Issue Date: 12/08/2009	Policy #: EHS-100.05								
X-Ray Machine Inspection Policy									

1.0 PURPOSE:

To ensure all X-ray machines at LSUHSC are properly inspected in accordance with radiation protection regulations and campus radiation safety committee guidance.

2.0 SCOPE:

This policy informs the X-Ray machine operator/inspectors/owners what items are to be addressed during the required inspection and before a Louisiana State Department of Environmental Quality (DEQ) X-ray machine audit.

3.0 RESPONSIBILITIES:

3.1 Radiation Safety Officer (RSO) shall:

- Perform all dental and veterinarian X-ray machine inspections every 3 years.
- Verify that all medical X-ray machine inspections are completed annually.
- Maintain all inventory records of X-Ray devices and inspection results.

3.2 Department X-Ray Owners shall:

- Notify RSO of any new X-ray machine purchase to ensure proper DEQ registration.
- Notify RSO of any X-ray machine removal.
- Notify RSO of any X-ray relocations.
- Notify RSO of any X-ray repairs and/or modifications.

4.0 IMPLEMENTATION REQUIREMENTS:

4.1 Inspection Requirements:

- Inspections shall be performed at determined intervals (See Section 3.1) and after any X-ray machine relocation, modification, or repair.
- The following items are required in the immediate equipment area of every X-ray unit:
 - Copy of DEQ Registration License.
 - Copy of DRC-3 form.
 - Technique chart which indicates what time duration (ms) and milliamp (ma). values are used for different patient sizes and/or body parts.
 - Warning Label posted on machine.



- Exposure indication of a visual or audible type when X-rays are produced.
- Shielding to protect the patient from scatter X-rays. A lead apron should be worn except for direct focused dental machines.
- Operator shall stand at least 12 feet from tube housing while making exposures or shall stand behind approved leaded glass shield.
- Ensure the DEQ registration license is correct (e.g., serial #, model #).
- Ensure the operator is located at least 12 feet from tube housing while making inspection exposures. Operator may stand behind glass leaded shield.
- X-Ray inspection specific tasks include recording the following:
 - 1. Exposure Duration (Time) Reproducibility test (use 4 timing tests)
 - 2. Exposure Reproducibility test (use 4 exposures made in one hour)
 - 3. Linearity test (if equipment allows choice of X-Ray current settings)
 - 4. Accuracy test (no more than 10% error from one reading to another)
 - 5. Use Appendix A, Form RS 04, X-ray Inspection Form, to record results.

5.0 **RECORDKEEPING:**

Copies of all inspections must be on file at the unit's location, with the RSO, and be monitored for the current fiscal year and the previous three fiscal years.

6.0 INSPECTIONS AND PROGRAM REVIEW:

This procedure shall be performed at the determined intervals (See Section 3.1) schedule or whenever there are any modifications or repairs to the X-ray machine.

7.0 **REFERENCE:**

LA DEQ Title 33, Part XV - Sections 603, 604 and 608

8.0 APPENDIX:

A. Radiation X-ray Inspection Form

LSU-HSC Radiation Safety Office X-Ray Machine Inspection

Facility Name			Facility Location				Date		
Building Name		Room # \$							
Unit Type(Dental, Medical, etc.)		State DEQ Registration #							
Manufacturer			Model				Ser	rial #	
	«Vp		meae				001	iai //	
		Г	Most	requer	t Exam Set	tina			
Person Interviewed			Most Frequent Exam Setting kVp mA			•	timo		
			κνρ		1117				
*Note:	Instruments used for measure	ments			Scatter S	Survey I	Instrument		
TIME	R/EXPOSURE REPRODUC	BILITY TEST	Bas	ed on se	ttings:	kVp)	mA	mSec
(T _{max}	− T _{min}) ≤0.1 Tavg (needs to be	e less than 10% err	or)	(E _{max}	(<i>−E</i> min)≤ 0.1	0 Eavo) (needs to	be less than 10)% error)
1	Dose Reading (mR)	Time (mS)		(kVp)			(Hvl)		
-									_
2									_
3									
									_
Avg									
	Southar Dadiation Massure	manta	_	Tabla 1	- values must be	arootor	than ana sha	wn bolow	1
	Scatter Radiation Measure	ements		Table 1	- values must be			All Other	1
Operator @ 12 feet= uR/ hr exp		OI	Design Derating Range	Measured Potential	Dental Intraoral Manufactured before 8/1/74 and		Diagnostic X-ray Half- Value Layer		
ι	uR per sec= x exp p	er sec		bystem	(kVp)		or before 12/1/80	(mm of Aluminum)	
=uR per exp		Below 51		30 40		N/A N/A	0.3	•	
				50		1.5	0.5		
		_		51		1.5	1.2		
		5	1 to 70	60 70		1.5 1.5	1.3	•	
				71	2.1		2.1	1	
				80		2.3	2.3		
				<u> </u>		2.5 2.7	2.5	4	
			Δ1	bove 70	100		3.0	3.0	-
				5000 70	120		3.2	3.2	1
					130		3.5	3.5	1
					140		3.8	3.8	
					150		4.1	4.1	J
<u>Facilit</u>	ty Equipment and Design								
_	Registration Certificate Adequate Signs Postec	I		RC 3 Po nielding/	sted Aprons prov	vided	Tech	nique Chart	Posted
	Dead man Type Expos	ure Switch	E×	posure	Switch Loc	ated A	Adequate <u>:</u>	> 10 feet fror	n tube
	NO DEVIATIONS FOU	ND		DEVIA	TIONS FO	UND			
<u>Comr</u>	nents:								
Inspe	ection performed by _				DATE				